

RELIANCE ERGONOMIC TABLETOP

Designed with your comfort in mind

Look closer. See further.

 HAAG-STREIT
RELIANCE




IMPROVED COMFORT

Ergonomic Tabletop

Key features & benefits:

- The innovatively designed table cut-out allows the doctor to have improved comfort and position at the slit lamp.
- Uniquely designed 20° armrests for enhanced comfort during and after a patient exam.
- Efficient cable management designed to keep cables and wires neatly organized and out of sight.
- Can be retrofitted to an existing 7900 or 7800 Reliance Instrument Stand to optimize the doctor's ergonomic positioning for improved comfort.



RELIANCE ERGONOMIC TABLETOP

Cut-out is approximately 2-3/4" deep. Front cut-out allows doctors to get closer to slit lamp for improved posture.

Dimension between the doctor and patient:

- 10-5/8" for BI/BM
- 12-5/8" for BP/BQ

20° bend with smooth transition for ergonomic arm placement.

Models	Dimensions (inches)	Weight (lb)
BQ/BP Rheostat	15.44 deep x 24 wide	17.90
BM/BI Rheostat	13.44 deep x 24 wide	17.20
BQ/BP Imaging	15.44 deep x 24 wide	18.15
BI Imaging	13.44 deep x 24 wide	17.45

Surface finish color, all models: White
Ergonomic Tabletop pairs with Haag-Streit LED Slit Lamps only.



HAAG-STREIT USA, Inc.

3535 Kings Mills Road
Mason, OH 45040
USA
Phone +15133367255
www.haag-streit-usa.com