

THE POWER OF TWO LEADERS **THE ADVANTAGE FOR YOUR PRACTICE**



When you partner with Haag-Streit USA and Essilor Instruments,
you align with two global leaders supporting a practice
built for today's clinical demands and future innovation.

Designed to Adapt

As your practice evolves - from traditional workflows to digital integration.

WORKPLACE SOLUTION

VERSATILE. ROBUST. BUILT TO LAST.

- Reliance-built exam chair & stand – engineered for daily clinical use and long-term durability
- Made in the USA by Haag-Streit Reliance
- Power hydraulic base with assisted recline, including full-flat positioning for patient comfort

Choose the lane that fits your needs:

- **Manual lane:**
Pair with a Haag-Streit slit lamp for proven optical performance and mechanical precision, backed by a 10-year warranty
- **Digital lane:**
Upgrade to Haag-Streit's Elara 900 slit lamp and Essilor Instruments' Vision-R 700 to introduce digital documentation, guided refraction, and a more efficient exam workflow



ELARA 900

SLIT LAMP IMAGING

See more. Work smarter. Feel better.

- A superior clinical view utilizing legendary Haag-Streit optics
- Customized workflow with efficient digital presets
- Optimized ergonomics for the examiner and patient
- Revolutionary projector light source, able to emit colors across the entire visible spectrum
- High-resolution 3D video and image capture in 4K



VISION-R 700

REFRACTION SYSTEM

- Exclusive optical module permitting multiple incremental power changes to secure rapid 0.25D proven accuracy
- Simultaneous lens changes for a more direct and controlled refraction method
- Smart programs for easy refractions
- Instant power changes for easy comparisons and patient confidence
- Improved patient experience thanks to a wider field of vision, silent lens changes, and an easy response method



Anterior & Posterior Imaging

CELLVIEW WRI-1

RETINAL IMAGER

High-quality ultra-widefield fundus imaging


- Captures a retinal image and far into the periphery, up to 133° in a single-capture, or an up to 200° auto-stitched image
- Provides clear and high-definition images through small pupils, most cataracts, and other media opacities
- Full visible spectrum LED array, with true color and infrared imaging
- Easy-to-operate for nominal staff training



CORNEA 550+

MULTIFUNCTIONAL CORNEAL TOPOGRAPHY

- Corneal topography
- Customized contact lens fitting
- Pupillography
- Advanced tear film analysis
- Meibography



LIMITED TIME SPECIAL PRICING FOR VISION ASSOCIATES

Contact:



Talia Bruno, Sr. Director of National Accounts
talia.bruno@essilorinstrumentsusa.com
c: 770.686.8937



Dyanna Chanthaboury, Key Accounts Director
dyanna.chanthaboury@haag-streit.com
c: 513.767.8759